

# M-Series Stretchers

**stryker**<sup>®</sup>



SM104 with Fifth Wheel



SM204 with Big Wheel



SM304 with Zoom<sup>®</sup> Drive

Medical

6300 South Sprinkle Road  
Kalamazoo, Mich. 49001 USA  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

## Standard Features

- Thermoformed ABS base with extra capacity storage tray
- 700 lbs weight capacity
- 26" wide patient surface with 3" (8 cm) Enhanced Comfort Mattress
- Less than 21" low height
- Integrated utility tray
- Litter mounted oxygen bottle holder
- Integrated transfer board
- Powerwashable
- Four-wheel steel-ring brake system with dual end activators
- Dual side-mounted foot control hydraulics with uni-lower pedal
- Glideaway™ siderails with integrated push/pull handles
- Four integrated IV receptacles
- Trend/reverse Trend.
- 8" (20 cm) Omni Surface™ casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual end brake & steer controls

SM 104: All Standard Features Plus

- Fifth Wheel steering

SM 204: All Standard Features Plus

- Big Wheel steering

SM 304: All Standard Features Plus

- Zoom Motorized Drive System

## Optional Features

- 30" width patient surface
- 4" (10 cm) Enhanced Comfort Mattress
- 4" (10 cm) or 5" (13 cm) Ultra Comfort Mattress
- Four-sided brake & steer controls (SM104 & SM204 only)
- Pop-up steering handles (SM104 & SM204 only)
- No-crank hydraulic knee gatch
- LiftAssist™ backrest
- Integrated scale system

## Specifications

	26" Litter		30" Litter	
• <b>Model Number</b>	<b>1005 (SM104)</b>		<b>1005 (SM104)</b>	
	<b>1015 (SM204)</b>		<b>1015 (SM204)</b>	
	<b>1025 (SM304)</b>		<b>1025 (SM304)</b>	
• <b>Overall Length</b>				
SM104 & SM204	84"	(213 cm)	84"	(213 cm)
SM304	85"	(215 cm)	85"	(215 cm)
• <b>Overall Width</b>				
Siderails Up	33.5"	(85 cm)	37"	(94 cm)
Siderails Down	30.25"	(77 cm)	30.75"	(78 cm)
• <b>Weight Capacity</b>	700 lbs	(317 kg)	700 lbs	(317 kg)
• <b>Height Range</b>				
High	34"	(86 cm)	34"	(86 cm)
Low (SM104 & SM204)	20.75"	(53 cm)	20.75"	(53 cm)
Low (SM304)	23"	(58 cm)	23"	(58 cm)
• <b>Litter Positioning</b>				
Backrest	0° - 90°		0° - 90°	
Knee Gatch	40°		40°	
Trend./reverse Trend. (SM104 & SM204)	±16°		±16°	
Trend./reverse Trend. (SM304)	±15°		±15°	
Maximum Angle	87.5°		87.5°	
• <b>Patient Surface</b>	26" x 74.25"	(66 x 189 cm)	30" x 74.25"	(76 x 189 cm)
• <b>Siderails</b>	14" x 57.75"	(36 x 147 cm)	14" x 57.75"	(36 x 147 cm)
• <b>Caster Diameter</b>	8"	(20 cm)	8"	(20 cm)

Stryker reserves the right to change specifications without notice.

## Accent Colors

- Red / Yellow / Light Blue / Purple / Gray / Dark Blue

## Warranty

- Two years parts, labor and travel